

TRAVEL TIPS FOR YOUR PUMP

- **When traveling with your pump it is recommended that you pack your pump in either a small suitcase or small duffel bag with some type of protective padding.**
- **We suggest using foam padding or bubble wrap, *which you can purchase from any local arts and craft store.***
- **It is imperative that the entire pump is completely wrapped so that it can withstand any chance of being bounced around.**
- **We advise that you place a visible tag or label stating “Handle with Care, Medical Equipment.”**
 - **We can provide you with some literature explaining what the pump is used for in the event that airport security has questions pertaining to the pump. Just give us a call a couple weeks before your trip and we will be happy to mail it right out to you.**
- **It is completely safe for any pump to go through the X-ray machines at the airport.**